

GST-tag, HRP-conjugated (mouse monoclonal, clone 3E2)

Cat PA00180

Size 100 µg

Host

Mouse

Clonality

Monoclonal

Confirmed reactivity

Immunogen

GST tagged protein from *Schistosoma japonicum*.

Host

Mouse

Clonality

Monoclonal

Purity

Total IgG1. Protein G purified in PBS pH 7.4, with 50% glycerol and 0.09% sodium azide.

Format

Liquid

Storage

Store at -20°C. Avoid repeated freeze-thaw cycles.

Application

ELISA (ELISA), Western blot (WB)

Recommended dilution

1 : 800 (ELISA), 1: 1000 (WB)

Expected | apparent MW

26 kDa

Additional information

For in vitro research use only

Description

The gst tag (glutathione S-transferase) is a tag that is added to the protein of interest as a fusion protein and cannot be purified and detected. The label is 220 amino acids long, about 26 kDa, which is quite large compared to Myc or FLAG labels.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY