

# Pentatricopeptide repeat-containing protein At1g14470 (PCMP-A4), Recombinant Protein

Cat *RP08436*

---

## Species

*Arabidopsis thaliana* (Mouse-ear cress)

## Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At1g14470 (PCMP-A4) , partial

## Product Gene Name

PCMP-A4 recombinant protein

## Product Synonym Gene Name

PCMP-A4

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

E Coli or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

61,742 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

## Protein Family

Pentatricopeptide repeat-containing protein

## NCBI Accession #

NP\_172899.1

## NCBI GI #

15223763

## NCBI GenBank Nucleotide #

NM\_101314.2

## NCBI GenID

838009

## NCBI Official Full Name

Pentatricopeptide repeat (PPR) superfamily protein

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

---

# Pentatricopeptide repeat-containing protein At1g14470 (PCMP-A4), Recombinant Protein

Cat *RP08436*

---

## NCBI Official Symbol

AT1G14470

## NCBI Official Synonym Symbols

F14L17.25; F14L17\_25

## NCBI Protein Information

Pentatricopeptide repeat (PPR) superfamily protein

## UniProt Gene Name

PCMP-A4

## UniProt Protein Name

Pentatricopeptide repeat-containing protein At1g14470

## UniProt Primary Accession #

Q9M9R6

## UniProt Related Accession #

Q9M9R6

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**