

DUF21 domain-containing protein At4g14230 (CBSDUF2), Recombinant Protein

Cat *RP08633*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* DUF21 domain-containing protein At4g14230 (CBSDUF2), partial

Product Gene Name

CBSDUF2 recombinant protein

Product Synonym Gene Name

CBSDUF2

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

53,492 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

DUF21 domain-containing protein

NCBI Accession

NP_193159.3

NCBI GI

42566780

NCBI GenBank Nucleotide

NM_117500.6

NCBI GeneID

827064

NCBI Official Full Name

CBS domain protein with a domain protein (DUF21)

NCBI Official Symbol

AT4G14230

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

DUF21 domain-containing protein At4g14230 (CBSDUF2), Recombinant Protein

Cat *RP08633*

NCBI Official Synonym Symbols

DL3155C; FCAALL.147

NCBI Protein Information

CBS domain protein with a domain protein (DUF21)

UniProt Gene Name

CBSDUF2

UniProt Protein Name

DUF21 domain-containing protein At4g14230

UniProt Synonym Protein Names

CBS domain-containing protein CBSDUF2

UniProt Primary Accession

Q4V3C7

UniProt Secondary Accession

O23281

UniProt Related Accession

Q4V3C7

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY