

Pisum sativum Lectin (PSA/PSL) - Separopore® 4B MicroPass™



Cat PL00045

Size 300 µL

Sugar Specificity

Description

MicroPass™ microcolumns are prepared by coupling Separopore® 4B with affinity-purified (AffyPure™) Pisum sativum lectin (PSA/PSL). These microcolumns are used for mass spectrometry and proteomics analysis of samples

Shelf Life

1 year

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY