

# Putative cysteine-rich repeat secretory protein 35 (CRRSP35), Recombinant Protein

Cat RP08645

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 1

mg (E-Coli)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)

## Species

Arabidopsis thaliana (Mouse-ear cress)

## Full Product Name

Recombinant Arabidopsis thaliana Putative cysteine-rich repeat secretory protein 35 (CRRSP35)

## Product Gene Name

CRRSP35 recombinant protein

## Product Synonym Gene Name

CRRSP35

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Sequence

LNLTNEYLNH KCFVSEGKYK HGDKYENNLN VLNKNVLSYD LTSGFLHVSH GEGPDSVTII LQCRGDSFGS  
NCRSCYTTAI DGFHRRRCQRN KGGIIFYDQC FLVISTIKPQ LPRKIDFKNT FSMHNPNNVS SEPGSFDKMT  
RDFLYKLVQK ASYPTVVEHQ STYYAAGEKK LGKRKLYAMM QCASDILQCK VCLEWCIREL PKCCDGKQGG  
RILGMSCNLR YELYPFLRR

## Sequence Positions

30-258, Full length protein

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

E Coli or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

29,634 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

## Protein Family

Putative cysteine-rich repeat secretory protein

## NCBI Accession #

NP\_188843.2

## NCBI GI #

42565101

## NCBI GenBank Nucleotide #

NM\_113101.3

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

# Putative cysteine-rich repeat secretory protein 35 (CRRSP35), Recombinant Protein

Cat      *RP08645*

Size      0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 1

mg (E-Coli)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)

## NCBI GeneID

821765

## NCBI Official Full Name

cysteine-rich repeat secretory protein, putative (DUF26)

## NCBI Official Symbol

AT3G22050

## NCBI Protein Information

cysteine-rich repeat secretory protein, putative (DUF26)

## UniProt Gene Name

CRRSP35

## UniProt Protein Name

Putative cysteine-rich repeat secretory protein 35

## UniProt Primary Accession #

Q9LRK2

## UniProt Related Accession #

Q9LRK2

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**