

Pentatricopeptide repeat-containing protein At3g05340 (PCMP-E83), Recombinant Protein

Cat *RP08747*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At3g05340 (PCMP-E83) , partial

Product Gene Name

PCMP-E83 recombinant protein

Product Synonym Gene Name

PCMP-E83

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

73,128 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_187185.2

NCBI GI

42563517

NCBI GenBank Nucleotide

NM_111407.2

NCBI GenID

819698

NCBI Official Full Name

Tetratricopeptide repeat (TPR)-like superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At3g05340 (PCMP-E83), Recombinant Protein

Cat *RP08747*

NCBI Official Symbol

AT3G05340

NCBI Official Synonym Symbols

T12H1.31; T12H1_31

NCBI Protein Information

Tetratricopeptide repeat (TPR)-like superfamily protein

UniProt Gene Name

PCMP-E83

UniProt Protein Name

Pentatricopeptide repeat-containing protein At3g05340

UniProt Primary Accession

Q9MA85

UniProt Related Accession

Q9MA85

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY