

Pentatricopeptide repeat-containing protein At4g39530 (PCMP-E52), Recombinant Protein

Cat *RP08750*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At4g39530 (PCMP-E52) , partial

Product Gene Name

PCMP-E52 recombinant protein

Product Synonym Gene Name

PCMP-E52

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

93,794 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_195663.1

NCBI GI

15235115

NCBI GenBank Nucleotide

NM_120113.2

NCBI GenID

830107

NCBI Official Full Name

Tetratricopeptide repeat (TPR)-like superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At4g39530 (PCMP-E52), Recombinant Protein

Cat *RP08750*

NCBI Official Symbol

AT4G39530

NCBI Official Synonym Symbols

F23K16.160; F23K16_160

NCBI Protein Information

Tetratricopeptide repeat (TPR)-like superfamily protein

UniProt Gene Name

PCMP-E52

UniProt Protein Name

Pentatricopeptide repeat-containing protein At4g39530

UniProt Primary Accession

Q9SVA5

UniProt Related Accession

Q9SVA5

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY