

Pentatricopeptide repeat-containing protein At2g25580 (PCMP-H75), Recombinant Protein

Cat *RP08752*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At2g25580 (PCMP-H75) , partial

Product Gene Name

PCMP-H75 recombinant protein

Product Synonym Gene Name

PCMP-H75

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

69,906 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_001324433.1

NCBI GI

1063704830

NCBI GenBank Nucleotide

NM_001336003.1

NCBI GenID

817097

NCBI Official Full Name

Tetratricopeptide repeat (TPR)-like superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At2g25580 (PCMP-H75), Recombinant Protein

Cat *RP08752*

NCBI Official Symbol

MEF8

NCBI Official Synonym Symbols

F3N11.3; F3N11_3; mitochondrial editing factor 8

NCBI Protein Information

Tetratricopeptide repeat (TPR)-like superfamily protein

UniProt Gene Name

PCMP-H75

UniProt Protein Name

Pentatricopeptide repeat-containing protein At2g25580

UniProt Primary Accession

Q680H3

UniProt Secondary Accession

Q9SLA5

UniProt Related Accession

Q680H3

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY