

Pentatricopeptide repeat-containing protein At1g04840 (PCMP-H64), Recombinant Protein

Cat *RP08757*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At1g04840 (PCMP-H64) , partial

Product Gene Name

PCMP-H64 recombinant protein

Product Synonym Gene Name

PCMP-H64

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

75,105 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_171976.1

NCBI GI

15220333

NCBI GenBank Nucleotide

NM_100362.2

NCBI GeneID

839402

NCBI Official Full Name

Tetratricopeptide repeat (TPR)-like superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At1g04840 (PCMP-H64), Recombinant Protein

Cat *RP08757*

NCBI Official Symbol

AT1G04840

NCBI Official Synonym Symbols

F13M7.17; F13M7_17

NCBI Protein Information

Tetratricopeptide repeat (TPR)-like superfamily protein

UniProt Gene Name

PCMP-H64

UniProt Protein Name

Pentatricopeptide repeat-containing protein At1g04840

UniProt Primary Accession

Q9MAT2

UniProt Related Accession

Q9MAT2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY