

Pentatricopeptide repeat-containing protein At5g50990 (PCMP-H59), Recombinant Protein

Cat *RP08760*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At5g50990 (PCMP-H59) , partial

Product Gene Name

PCMP-H59 recombinant protein

Product Synonym Gene Name

PCMP-H59

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

60,036 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_199912.2

NCBI GI

145359106

NCBI GenBank Nucleotide

NM_124478.3

NCBI GenID

835172

NCBI Official Full Name

Tetratricopeptide repeat (TPR)-like superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At5g50990 (PCMP-H59), Recombinant Protein

Cat *RP08760*

NCBI Official Symbol

AT5G50990

NCBI Official Synonym Symbols

K3K7.15; K3K7_15

NCBI Protein Information

Tetratricopeptide repeat (TPR)-like superfamily protein

UniProt Gene Name

PCMP-H59

UniProt Protein Name

Pentatricopeptide repeat-containing protein At5g50990

UniProt Primary Accession

Q9FI49

UniProt Related Accession

Q9FI49

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY