

Pentatricopeptide repeat-containing protein At4g39952, mitochondrial (PCMP-E98), Recombinant Protein

Cat *RP08761*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At4g39952, mitochondrial (PCMP-E98) , partial

Product Gene Name

PCMP-E98 recombinant protein

Product Synonym Gene Name

PCMP-E98

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

85,828 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_680777.2

NCBI GI

334187347

NCBI GenBank Nucleotide

NM_148411.3

NCBI GeneID

830155

NCBI Official Full Name

Pentatricopeptide repeat (PPR) superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At4g39952, mitochondrial (PCMP-E98), Recombinant Protein

Cat *RP08761*

NCBI Official Symbol

AT4G39952

NCBI Protein Information

Pentatricopeptide repeat (PPR) superfamily protein

UniProt Gene Name

PCMP-E98

UniProt Protein Name

Pentatricopeptide repeat-containing protein At4g39952, mitochondrial

UniProt Primary Accession

Q3E9N1

UniProt Related Accession

Q3E9N1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY