

Pentatricopeptide repeat-containing protein At5g65570 (PCMP-H47), Recombinant Protein

Cat *RP08764*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At5g65570 (PCMP-H47) , partial

Product Gene Name

PCMP-H47 recombinant protein

Product Synonym Gene Name

PCMP-H47

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

83,681 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_201360.1

NCBI GI

15239085

NCBI GenBank Nucleotide

NM_125955.2

NCBI GenID

836683

NCBI Official Full Name

Tetratricopeptide repeat (TPR)-like superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At5g65570 (PCMP-H47), Recombinant Protein

Cat *RP08764*

NCBI Official Symbol

AT5G65570

NCBI Official Synonym Symbols

K21L13.8; K21L13_8

NCBI Protein Information

Tetratricopeptide repeat (TPR)-like superfamily protein

UniProt Gene Name

PCMP-H47

UniProt Protein Name

Pentatricopeptide repeat-containing protein At5g65570

UniProt Primary Accession

Q9LSL8

UniProt Related Accession

Q9LSL8

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY