

Pentatricopeptide repeat-containing protein At3g14730 (PCMP-E31), Recombinant Protein

Cat *RP08800*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At3g14730 (PCMP-E31) , partial

Product Gene Name

PCMP-E31 recombinant protein

Product Synonym Gene Name

PCMP-E31

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

72,711 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_188091.1

NCBI GI

15231917

NCBI GenBank Nucleotide

NM_112334.2

NCBI GenID

820701

NCBI Official Full Name

Pentatricopeptide repeat (PPR) superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At3g14730 (PCMP-E31), Recombinant Protein

Cat *RP08800*

NCBI Official Symbol

AT3G14730

NCBI Protein Information

Pentatricopeptide repeat (PPR) superfamily protein

UniProt Gene Name

PCMP-E31

UniProt Protein Name

Pentatricopeptide repeat-containing protein At3g14730

UniProt Primary Accession

Q9LUC2

UniProt Related Accession

Q9LUC2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY