

F-box protein At3g20690 (At3g20690), Recombinant Protein

Cat RP09326

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Species Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant Arabidopsis thaliana F-box protein At3g20690 (At3g20690)

Product Gene Name

At3g20690 recombinant protein

Product Synonym Gene Name

At3g20690

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MMMSDLPHDL VEEILSRLPL ISLKAMRSTC KTNVLSKHR SFANKHIGNV TASGKRDLIM IKDCKVYSIG
VNLHGIQNNN NIIDLSIKNK GILHLENLDL IFKEIFNVFH CNGLLLLYGS ITDDSIRLVV CNPYWGKREW
VKRINNFAF DKLAFGYDKS CGCHKILRVF DSYPNRLEIY NLRSNSRMVS SFPLKWDIAF MQDGVSLKGN
TYWYAKDRRS EDCYLICFDF TRERFGPRLP LPQPFIDEFL GSLYVVREEK LAVLVKRWRG GKSEIDIWVT
NKIEPDEVSX SKFLKVDNTP RFFLTEGYLI YDEEKVVVFF NKEETEGLS TVHIAYIVGE NGYFRRADLR
EPLYQPLCRL VCSYVPSSVK

Sequence Positions

1-370, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

13,125 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

F-box protein

NCBI Accession

NP_188705.1

NCBI GI

15232336

NCBI GenBank Nucleotide

NM_112960.1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

F-box protein At3g20690 (At3g20690), Recombinant Protein

Cat *RP09326*

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-
Cell)

NCBI GeneID

821616

NCBI Official Full Name

F-box and associated interaction domains-containing protein

NCBI Official Symbol

AT3G20690

NCBI Protein Information

F-box and associated interaction domains-containing protein

UniProt Gene Name

At3g20690

UniProt Protein Name

F-box protein At3g20690

UniProt Primary Accession

Q1PEN8

UniProt Secondary Accession

Q9LHQ3

UniProt Related Accession

Q1PEN8

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY