

B3 domain-containing protein At5g57720 (At5g57720), Recombinant Protein

Cat RP09811

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-
Arabidopsis thaliana (Mouse-ear cress)
Cell)

Full Product Name

Recombinant Arabidopsis thaliana B3 domain-containing protein At5g57720 (At5g57720)

Product Gene Name

At5g57720 recombinant protein

Product Synonym Gene Name

At5g57720

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MVLENGFPPY PDFLKIFNSH EDSQLLVIPR SYNRYYPNPL PQTAVLKNPE GRFWNVQWTK SQEVIISLQE
GWVKFVKDNG LIDRDFLLFT YDGSRSFWVR IHRNGLPLEP TAPIKIQEIS DDEDETNGDG DPHMEEEGDT
DENMIVSLSL GSSDEVGDDD DDDYDTAICD EVNKASGSSK KGRVTRKHRD DSLASITPQI FLADPNNPFF
ISTSSCSRRI LVIARQVIKD YGLNFDGTVN LIDGFGELTR KVGKWKDRVV IYKWNEMFTR NKVKQGDVII
CEIIREEDVV RSIKVHFVKN

Sequence Positions

1-300, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

34,399 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week.
Repeated freezing and thawing is not recommended.

Protein Family

B3 domain-containing protein

NCBI Accession

NP_200580.1

NCBI GI

15242851

NCBI GenBank Nucleotide

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

B3 domain-containing protein At5g57720 (At5g57720), Recombinant Protein

Cat RP09811

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-
Cell)
NCBI GeneID
835879

NCBI Official Full Name

AP2/B3-like transcriptional factor family protein

NCBI Official Symbol

AT5G57720

NCBI Official Synonym Symbols

MRI1.8; MRI1_8

NCBI Protein Information

AP2/B3-like transcriptional factor family protein

UniProt Gene Name

At5g57720

UniProt Protein Name

B3 domain-containing protein At5g57720

UniProt Primary Accession

Q9FHH1

UniProt Related Accession

Q9FHH1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY