

Pentatricopeptide repeat-containing protein At1g15480, mitochondrial (At1g15480), Recombinant Protein

Cat *RP09817*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g15480, mitochondrial (At1g15480) , partial

Product Gene Name

At1g15480 recombinant protein

Product Synonym Gene Name

At1g15480

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

67,532 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Pentatricopeptide repeat-containing protein

NCBI Accession

NP_173001.2

NCBI GI

145335742

NCBI GenBank Nucleotide

NM_101417.4

NCBI GeneID

838118

NCBI Official Full Name

Tetratricopeptide repeat (TPR)-like superfamily protein

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Pentatricopeptide repeat-containing protein At1g15480, mitochondrial (At1g15480), Recombinant Protein

Cat *RP09817*

NCBI Official Symbol

AT1G15480

NCBI Official Synonym Symbols

F9L1.43; F9L1_43

NCBI Protein Information

Tetratricopeptide repeat (TPR)-like superfamily protein

UniProt Gene Name

At1g15480

UniProt Protein Name

Pentatricopeptide repeat-containing protein At1g15480, mitochondrial

UniProt Primary Accession

Q9XI21

UniProt Related Accession

Q9XI21

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
