

# Pentatricopeptide repeat-containing protein At1g07590, mitochondrial (At1g07590), Recombinant Protein

Cat *RP09830*

---

## Species

*Arabidopsis thaliana* (Mouse-ear cress)

## Full Product Name

Recombinant *Arabidopsis thaliana* Pentatricopeptide repeat-containing protein At1g07590, mitochondrial (At1g07590) , partial

## Product Gene Name

At1g07590 recombinant protein

## Product Synonym Gene Name

At1g07590

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

*E. coli* or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

60,849 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

## Protein Family

Pentatricopeptide repeat-containing protein

## NCBI Accession #

NP\_172238.1

## NCBI GI #

15222438

## NCBI GenBank Nucleotide #

NM\_100632.3

## NCBI GeneID

837272

## NCBI Official Full Name

Tetratricopeptide repeat (TPR)-like superfamily protein

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Pentatricopeptide repeat-containing protein At1g07590, mitochondrial (At1g07590), Recombinant Protein

Cat *RP09830*

---

## NCBI Official Symbol

AT1G07590

## NCBI Official Synonym Symbols

F22G5.3; F22G5\_3

## NCBI Protein Information

Tetratricopeptide repeat (TPR)-like superfamily protein

## UniProt Gene Name

At1g07590

## UniProt Protein Name

Pentatricopeptide repeat-containing protein At1g07590, mitochondrial

## UniProt Primary Accession #

Q940Q2

## UniProt Secondary Accession #

Q9LNY2

## UniProt Related Accession #

Q940Q2

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**