

F-box protein At1g67130 (At1g67130), Recombinant Protein

Cat RP09904

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 0.1 mg (Mammalian-Cell)/ 1
mg (Yeast)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Species Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant Arabidopsis thaliana F-box protein At1g67130 (At1g67130)

Product Gene Name

At1g67130 recombinant protein

Product Synonym Gene Name

At1g67130

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MNRGETLDSI PTDLILDILS RLPTKSIARF HCVSKLWSSM LASQDFTRLF VNRSSSNPRL LLGIVRGGEW
SFYSSPQPKN PYEKSSLEVA ADFHMKLSEI KQHQQHGTSY AFGLIYLRTV WMSKEGNYEV CVICKPSTGQ
YYAILPPQPS GIFGDLGFDP IGKQFKVLVL NILDNHYILT LGTENDKWRS IQSSLRYRPC GQSPICINGV
LYYIAYDTQD SSNDVIGCFD VRFEKFKFFH VNHDMMVKCFF ELINYKDKLG GFPLELSMWV LEDLEKEEWS
KYAYTLKPDN NVVKVNYNLS FVGVTTRGEI VLAKMYACKP FYVFYYPEN NTLLSVEIQG VGEDSEWFSN
YQRVYVFVEH VKDLHQFDNT KTSINLPEQK RKPTSISISS KYDDHVRTTL ISSRKNQQVT TVSRPQQDRR
TTNKFSALCL LVMKNSQVSK PDCFSWSSAT MLYL

Sequence Positions

1-454, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

20,702 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week.
Repeated freezing and thawing is not recommended.

Protein Family

F-box protein

NCBI Accession

NP_176884.1

NCBI GI

15219833

NCBI GenBank Nucleotide

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

F-box protein At1g67130 (At1g67130), Recombinant Protein

Cat *RP09904*

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 0.1 mg (Mammalian-Cell)/ 1

mg (Yeast)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)

NCBI GeneID

843033

NCBI Official Full Name

F-box family protein

NCBI Official Symbol

AT1G67130

NCBI Official Synonym Symbols

F5A8.14; F5A8_14

NCBI Protein Information

F-box family protein

UniProt Gene Name

At1g67130

UniProt Protein Name

F-box protein At1g67130

UniProt Primary Accession

Q3ECH0

UniProt Secondary Accession

F4HRR9

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY