

Putative F-box protein At1g67390 (At1g67390), Recombinant Protein

Cat RP09974

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Yeast)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Species Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant Arabidopsis thaliana Putative F-box protein At1g67390 (At1g67390)

Product Gene Name

At1g67390 recombinant protein

Product Synonym Gene Name

At1g67390

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MSQFPLSILK TSQPHHQPKS MNAPRSSDAR GKRILIDPID DRISKLPDDV LVMILASLST EDALKTSVLS
TRWKNVWKQV PYLHFDLLSA TYESGLIAAL SNSVAESITQ VINNHNGHLM GCSIHHFAHH CVNGVLENWI
RLLTLEKNTR FLSLCNKLGK GRRTNVLRF S PNTFSHPKLA TLFLRRYDLV TAHAFKECEN LKILKLEGIL
AEVDVFNTVI ASCPSLKVLV LNCMWYNEKT CLKIHNNNLK VLHLDCPHVD CIDVSAALLD IFSIFYFAFA
KEQNFVIKAP RLLFNNWTKD LEKVPNMNYN ITSHAQEKKN IGHEFVVSGD ASYLQKLKSL KVVVDVTNSR
EVHMLRDILV AWNGVLEELH ILFKDNNVPN SKEDSESSIG GTQKKKWEEA NLFPNADFRV EVMWMFDFNG
SNKKQFALAS RFVTQGTVMK KMMIKTSSFF ANEKLDNEAV VAKLKELPKG NENLSIECF

Sequence Positions

1-479, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

54,474 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Putative F-box protein

NCBI Accession

NP_176906.1

NCBI GI

15220428

NCBI GenBank Nucleotide

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Putative F-box protein At1g67390 (At1g67390), Recombinant Protein

Cat *RP09974*

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 0.1 mg (Mammalian-Cell)/ 1

mg (Yeast)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)

NCBI GeneID

843059

NCBI Official Full Name

F-box family protein

NCBI Official Symbol

AT1G67390

NCBI Official Synonym Symbols

F1N21.21

NCBI Protein Information

F-box family protein

UniProt Gene Name

At1g67390

UniProt Protein Name

Putative F-box protein At1g67390

UniProt Primary Accession

Q9FYF4

UniProt Related Accession

Q9FYF4

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY