

Cell number regulator 5 (CNR5), Recombinant Protein

Cat *RP10506*

Species

Zea mays (Maize)

Full Product Name

Recombinant Zea mays Cell number regulator 5 (CNR5), partial

Product Gene Name

CNR5 recombinant protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

19,738 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Cell number regulator

NCBI Accession

NP_001132085.1

NCBI GI

212723940

NCBI GenBank Nucleotide

NM_001138613.1

NCBI GeneID

100193499

NCBI Official Full Name

cell number regulator 5

NCBI Official Symbol

CNR05

NCBI Official Synonym Symbols

CNR5; ZmCNR05; GRMZM5G892035

NCBI Protein Information

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Cell number regulator 5 (CNR5), Recombinant Protein

Cat *RP10506*

cell number regulator 5

UniProt Gene Name

CNR5

UniProt Protein Name

Cell number regulator 5

UniProt Synonym Protein Names

ZmCNR05

UniProt Primary Accession

B4FF80

UniProt Secondary Accession

B6TTK2

UniProt Related Accession

B4FF80

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY