

Cell number regulator 10 (CNR10), Recombinant Protein

Cat *RP10582*

Species

Zea mays (Maize)

Full Product Name

Recombinant Zea mays Cell number regulator 10 (CNR10), partial

Product Gene Name

CNR10 recombinant protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

16,660 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Cell number regulator

NCBI Accession

NP_001151128.1

NCBI GI

226508610

NCBI GenBank Nucleotide

NM_001157656.1

NCBI GeneID

100284761

NCBI Official Full Name

cell number regulator 10

NCBI Official Symbol

CNR10

NCBI Official Synonym Symbols

ZmCNR10; GRMZM2G325477

NCBI Protein Information

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Cell number regulator 10 (CNR10), Recombinant Protein

Cat *RP10582*

cell number regulator 10

UniProt Gene Name

CNR10

UniProt Protein Name

Cell number regulator 10

UniProt Synonym Protein Names

ZmCNR10

UniProt Primary Accession

D9HP26

UniProt Secondary Accession

B6TX86

UniProt Related Accession

D9HP26

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY