

# Formin-like protein 3 (FH3), Recombinant Protein

Cat *RP10919*

---

## Species

*Oryza sativa* subsp. *japonica* (Rice)

## Full Product Name

Recombinant *Oryza sativa* subsp. *japonica* Formin-like protein 3 (FH3), partial

## Product Gene Name

FH3 recombinant protein

## Product Synonym Gene Name

FH3

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

E Coli or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

134,361 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

## Protein Family

Defensin

## NCBI Accession #

Q7G6K7.2

## NCBI GI #

172046164

## NCBI Official Full Name

Formin-like protein 3

## UniProt Gene Name

FH3

## UniProt Protein Name

Formin-like protein 3

## UniProt Synonym Protein Names

OsFH3

## UniProt Primary Accession #

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

---

# Formin-like protein 3 (FH3), Recombinant Protein

Cat *RP10919*

---

Q7G6K7

## UniProt Secondary Accession #

Q0IZ88; Q7XH56; Q8W5K6

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**