

CASP-like protein Osl_39071, Recombinant Protein

Cat *RP11134*

Size *0.02 mg/ 0.1 mg/ 5x0.1 mg*

Species

Oryza sativa subsp. *indica* (Rice)

Full Product Name

Recombinant *Oryza sativa* subsp. *indica* CASP-like protein Osl_39071 (Osl_39071)

Product Gene Name

Osl_39071 recombinant protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MDLEKGGKKPS EQAAACRIMQ VKDKLITLQP VVRACVFLAT AVAAVIMGLN KQSYTTVVAI VGTRPVTQTF
TAKFKDTPAF VFFVIANAIA SGYNLMVLVT RRILQRRAS LSVHLLDMVI LTLLATGSAT AASMAQLGKN
GNLHARWNPI CDKFGSFCNH GGIALVSSFI GVALMLALNL LSAAANSPRS NVTGQ

Sequence Positions

1-195aa; full length protein

Format

Liquid containing glycerol

Host

Cell Free Expression

Molecular Weight

20,765 Da

Storage

Store at -20°C, for extended storage, conserve at -20°C or -80°C. Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

NCBI Accession

B8BM7.1

NCBI GI

341958535

NCBI Official Full Name

CASP-like protein 1B2

UniProt Gene Name

Osl_39071

UniProt Synonym Gene Names

OsCASPL1B2

UniProt Protein Name

CASP-like protein 1B2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

CASP-like protein Osl_39071, Recombinant Protein

Cat *RP11134*

Size *0.02 mg/ 0.1 mg/ 5x0.1 mg*

UniProt Entry Name

CSPL6_ORYSI

UniProt Primary Accession

B8BMY7

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
