

Zinc finger A20 and AN1 domain-containing stress-associated protein 11 (SAP11), Recombinant Protein

Cat *RP12301*

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 1

mg (E-Coli)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Species *Oryza sativa* subsp. japonica (Rice)

Full Product Name

Recombinant *Oryza sativa* subsp. japonica Zinc finger A20 and AN1 domain-containing stress-associated protein 11 (SAP11)

Product Gene Name

SAP11 recombinant protein

Product Synonym Gene Name

SAP11

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MAQREKKVEE PTELRAPEMT LCANSCGFPG NPATNNLCQN CFLAASASSS SSSAAASPST TSLPVFPVVE
KPRQAVQSSA AAVALVVER PTAGPVESSS KASRSSSVNR CHSCRRRVGL TGFRRCRCGEL YCGAHRYSDR
HDCSFDYKSA ARDAIAREN P VVRAAKIVRF

Sequence Positions

1-170, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

18,144 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Zinc finger A20 and AN1 domain-containing stress-associated protein

NCBI Accession

XP_015651039.1

NCBI GI

1002293396

NCBI GenBank Nucleotide

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Zinc finger A20 and AN1 domain-containing stress-associated protein 11 (SAP11), Recombinant Protein

Cat *RP12301*

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 1

mg (E-Coli)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
XM_015795593.1
NCBI GeneID
4345973

NCBI Official Full Name

zinc finger A20 and AN1 domain-containing stress-associated protein 11

NCBI Official Symbol

LOC4345973

NCBI Official Synonym Symbols

SAP11; ZFP60; OsSAP11; OJ1134_H03.11

NCBI Protein Information

zinc finger A20 and AN1 domain-containing stress-associated protein 11

UniProt Gene Name

SAP11

UniProt Synonym Gene Names

ZFP60; OsSAP11

UniProt Protein Name

Zinc finger A20 and AN1 domain-containing stress-associated protein 11

UniProt Primary Accession

Q84PD8

UniProt Secondary Accession

A3BUM6

UniProt Related Accession

Q84PD8

UniProt Comments

May be involved in environmental stress response.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY