

SRL2, Recombinant Protein

Cat *RP12395*

Size *0.05 mg/ 0.1 mg/ 5x0.1 mg*

Species

Oryza sativa subsp. *japonica* (Rice)

Full Product Name

Recombinant Rice SRL2 protein, N-His tag

Product Synonym Names

SRL2; Protein SEMI-ROLLED LEAF 2

Product Gene Name

SRL2 recombinant protein

Purity

>90% as determined by SDS-PAGE

Format

Lyophilized. Lyophilized from PBS pH 7.4, 0.2% NLS, 5% Trehalose, 5% Mannose, 1mM EDTA.

Host

E Coli

Molecular Weight

Predicted Molecular Mass: 114.6kDa

Storage

Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. Use a manual defrost freezer and avoid repeated freeze thaw cycles. In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Protein Family

Protein

UniProt Primary Accession

Q10MI0

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY