

4-coumarate--CoA ligase-like 4 (4CLL4), Recombinant Protein

Cat RP12622

Species

Oryza sativa subsp. japonica (Rice)

Full Product Name

Recombinant Oryza sativa subsp. japonica 4-coumarate--CoA ligase-like 4 (4CLL4), partial

Product Gene Name

4CLL4 recombinant protein

Product Synonym Gene Name

4CLL4

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

59,865 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

NCBI Accession

XP_015629747.1

NCBI GI

1002251034

NCBI GenBank Nucleotide

XM_015774261.1

NCBI GeneID

4331515

NCBI Official Full Name

4-coumarate--CoA ligase-like 4

NCBI Official Symbol

LOC4331515

NCBI Official Synonym Symbols

4CLL4

NCBI Protein Information

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

4-coumarate--CoA ligase-like 4 (4CLL4), Recombinant Protein

Cat *RP12622*

4-coumarate--CoA ligase-like 4

UniProt Gene Name

4CLL4

UniProt Protein Name

4-coumarate--CoA ligase-like 4

UniProt Primary Accession

Q10S72

UniProt Secondary Accession

B7EYA8

UniProt Related Accession

Q10S72

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY