

Arachis hypogaea Lectin (PNA) - Separopore® 4B

Cat PL00094

Size 10 mL/25 mL/50 mL/2 mL/5 mL

Sugar Specificity

Gal β 3GalNAc

Description

Affinity purified peanut lectin (PNA) is fixed on Separopore®4B by CNBr coupling. The final product is stabilized in the buffer to maintain maximum activity. The gel can be filled with 1 ml or 5 ml.

Abbreviation

PNA

Material Source

Peanut

Conjugate

Separopore® 4B

Purity

High Purity Grade

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

Gal β 3GalNAc

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Molecular Weight

N/A

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY