

Putative cyclin-F1-2 (CYCF1-2), Recombinant Protein

Cat RP13665

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Species Oryza sativa subsp. japonica (Rice)

Full Product Name

Recombinant Oryza sativa subsp. japonica Putative cyclin-F1-2 (CYCF1-2)

Product Gene Name

CYCF1-2 recombinant protein

Product Synonym Gene Name

CYCF1-2

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MFPGYSGGGG VPLPQMDADT YVRTIAAMPP HPLAPPPDSP RTPHTYVGFL PVFGDLPPLT GAVLQEPVPV
PPEQRADQPV AVATENSAPT RPQLCAPYDD EVEATLRAME TNPAERPSPY FLETTQGGRM SALVRASMIA
FMGEFSRKNK LADGTLQRAA YFLDRYLSVT PESDDALQLR LVGATAVFLA AKYEDQYTLR KIDASMVAAR
CGYTSETRHK MVSIMETEML AALGFNLGGP TAYTFVEHFT RYYGDGEEEEK QLKEAAHRVA DGTLLTYGFH
RYLPSMVAAS SIFLARLHEL GHEPWSNDLA ELTGKAILD MGCVCDIYSQ IACPRFALLQ GKNKRVCKST
SMQGMNMMMA

Sequence Positions

1-359, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

39,379 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Putative cyclin

NCBI Accession

Q6K8S5.1

NCBI GI

75291089

NCBI Official Full Name

Putative cyclin-F1-2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Putative cyclin-F1-2 (CYCF1-2), Recombinant Protein

Cat *RP13665*

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-
CYCF1-2 Cell)

UniProt Gene Name

UniProt Synonym Gene Names

CycF1;2

UniProt Protein Name

Putative cyclin-F1-2

UniProt Primary Accession

Q6K8S5

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY