

Peroxisomal membrane protein 11-3 (PEX11-3), Recombinant Protein

Cat *RP13738*

Species

Oryza sativa subsp. *japonica* (Rice)

Full Product Name

Recombinant *Oryza sativa* subsp. *japonica* Peroxisomal membrane protein 11-3 (PEX11-3), partial

Product Gene Name

PEX11-3 recombinant protein

Product Synonym Gene Name

PEX11-3

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

25,931 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Peroxisomal membrane protein

NCBI Accession

XP_015630855.1

NCBI GI

1002253235

NCBI GenBank Nucleotide

XM_015775369.1

NCBI GeneID

4332579

NCBI Official Full Name

peroxisomal membrane protein 11-3

NCBI Official Symbol

LOC4332579

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Peroxisomal membrane protein 11-3 (PEX11-3), Recombinant Protein

Cat *RP13738*

NCBI Official Synonym Symbols

PEX11-3; OsPEX11-3; OsJ_010110

NCBI Protein Information

peroxisomal membrane protein 11-3

UniProt Gene Name

PEX11-3

UniProt Protein Name

Peroxisomal membrane protein 11-3

UniProt Synonym Protein Names

OsPEX11-3; Peroxin-11-3

UniProt Primary Accession

Q10MN2

UniProt Secondary Accession

A0A0P0VWF8

UniProt Related Accession

Q10MN2

UniProt Comments

Involved in peroxisomal proliferation.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
