

Cyclin-T1-4 (CYCT1-1), Recombinant Protein

Cat *RP13756*

Species

Oryza sativa subsp. *japonica* (Rice)

Full Product Name

Recombinant *Oryza sativa* subsp. *japonica* Cyclin-T1-4 (CYCT1-1), partial

Product Gene Name

CYCT1-1 recombinant protein

Product Synonym Gene Name

CYCT1-1

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

61,391 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Cyclin

NCBI Accession

XP_015619900.1

NCBI GI

1002313032

NCBI GenBank Nucleotide

XM_015764414.1

NCBI GeneID

4352229

NCBI Official Full Name

cyclin-T1-4

NCBI Official Symbol

LOC4352229

NCBI Official Synonym Symbols

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Cyclin-T1-4 (CYCT1-1), Recombinant Protein

Cat *RP13756*

CYCT1-1

NCBI Protein Information

cyclin-T1-4

UniProt Gene Name

CYCT1-1

UniProt Synonym Gene Names

CycT1;4

UniProt Protein Name

Cyclin-T1-4

UniProt Primary Accession

Q2QQS5

UniProt Secondary Accession

B7EPW6

UniProt Related Accession

Q2QQS5

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY