

REF/SRPP-like protein Osl_017815 (Osl_017815), Recombinant Protein

Cat RP13856

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Species Oryza sativa subsp. indica (Rice)

Full Product Name

Recombinant Oryza sativa subsp. indica REF/SRPP-like protein Osl_017815 (Osl_017815)

Product Gene Name

Osl_017815 recombinant protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MADSGSDAPI SNRPEEEVTV EKTPEMEAAA EEERLRYLEF VQQAAAQVLV LAAAAYAYAK QGAGPLRPGV
DHVEGTVKAV VGPVYDRFHG VPLDLLKFLD RKGESVQEL DRRVPPVVKE APGLARSAAA EVRQAGLVGT
ATGLAKSAIA RAEPRARDLY TRYEPVAERK AAEAWAALNR LPLVPSVTRA VLPAAASLSA RYNTAVADGA
KRGSAVATYL PLVPTERLSR VFGYPLADAA ASPAPEMQPI PSQ

Sequence Positions

1-253, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

26,830 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

REF/SRPP-like protein

NCBI Accession

A2Y0H2.1

NCBI GI

187658937

NCBI Official Full Name

REF/SRPP-like protein Osl_017815

UniProt Gene Name

Osl_017815

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

REF/SRPP-like protein Osl_017815 (Osl_017815), Recombinant Protein

Cat *RP13856*

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-
Cell)

UniProt Primary Accession #

A2Y0H2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY