

Osl_21389, Recombinant Protein

Cat *RP13859*

Size 0.05 mg/ 0.1 mg/ 5x0.1 mg

Species

Oryza sativa subsp. indica (Rice)

Full Product Name

Recombinant Rice Osl_21389 protein, N-His Tag

Product Synonym Names

DsbD_2 domain-containing protein Osl_21389

Product Gene Name

Osl_21389 recombinant protein

Purity

>90% as determined by SDS-PAGE

Format

Lyophilized. Lyophilized from PBS pH 7.4, 0.2% NLS, 5% Trehalose, 5% Mannose, 1mM EDTA.

Host

E Coli

Molecular Weight

Predicted Molecular Mass: 22.4kDa

Storage

Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. Use a manual defrost freezer and avoid repeated freeze thaw cycles. In general, proteins are provided as lyophilized powder or frozen liquid. They are shipped out with dry ice or blue ice unless customers require otherwise.

UniProt Primary Accession

A2Y8K8

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY