

CASP-like protein Osl_03581 (Osl_03581), Recombinant Protein

Cat *RP13866*

Species

Oryza sativa subsp. indica (Rice)

Full Product Name

Recombinant Oryza sativa subsp. indica CASP-like protein Osl_03581 (Osl_03581), partial

Product Gene Name

Osl_03581 recombinant protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

21,771 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

NCBI Accession

B8A927.1

NCBI GI

341958557

NCBI Official Full Name

CASP-like protein 3A1

UniProt Gene Name

Osl_03581

UniProt Synonym Gene Names

OsCASPL3A1

UniProt Protein Name

CASP-like protein 3A1

UniProt Primary Accession

B8A927

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY