

CASP-like protein Os07g0442900 (Os07g0442900, LOC_Os07g26110), Recombinant Protein

Cat *RP14692*

Species

Oryza sativa subsp. *japonica* (Rice)

Full Product Name

Recombinant *Oryza sativa* subsp. *japonica* CASP-like protein Os07g0442900 (Os07g0442900, LOC_Os07g26110), partial

Product Gene Name

Os07g0442900 recombinant protein

Product Synonym Gene Name

Os07g0442900; LOC_Os07g26110

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

21,643 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

NCBI Accession

XP_015644614.1

NCBI GI

1002280623

NCBI GenBank Nucleotide

XM_015789128.1

NCBI GeneID

4343098

NCBI Official Full Name

CASP-like protein 1D1

NCBI Official Symbol

LOC4343098

NCBI Official Synonym Symbols

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

CASP-like protein Os07g0442900 (Os07g0442900, LOC_Os07g26110), Recombinant Protein

Cat *RP14692*

OsCASPL1D1; OsJ_023119

NCBI Protein Information

CASP-like protein 1D1

UniProt Gene Name

Os07g0442900

UniProt Synonym Gene Names

OsCASPL1D1

UniProt Protein Name

CASP-like protein 1D1

UniProt Primary Accession

Q6YT98

UniProt Secondary Accession

A0A0P0X5W9

UniProt Related Accession

Q6YT98

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY