

# B3 domain-containing protein LOC\_Os02g10420 (LOC\_Os02g10420), Recombinant Protein

Cat *RP14755*

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg  
(Yeast)/ 0.02 mg (Baculovirus)/ 1 mg (E-Coli)/ 0.02 mg

Species (Mammalian-Cell)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg  
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-  
Oryza sativa subsp. japonica (Rice)  
(Cell)

## Full Product Name

Recombinant Oryza sativa subsp. japonica B3 domain-containing protein LOC\_Os02g10420 (LOC\_Os02g10420)

## Product Gene Name

LOC\_Os02g10420 recombinant protein

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Sequence

MSAMLNENVP VEFQEAHGYA AREKVVLRLMR GRSWTVRLKH TKGRRPRRER AVLRYGWHRF  
CADNGLAVGD TCFFRALRSA GSGAGDVDDG DGDHVLSTV HKADGGDPLE

## Sequence Positions

1-110, Full length protein

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

E Coli or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

12,222 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week.  
Repeated freezing and thawing is not recommended.

## Protein Family

B3 domain-containing protein

## NCBI Accession #

XP\_015624348.1

## NCBI GI #

1002240183

## NCBI GenBank Nucleotide #

XM\_015768862.1

## NCBI GeneID

107276623

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# B3 domain-containing protein LOC\_Os02g10420 (LOC\_Os02g10420), Recombinant Protein

Cat *RP14755*

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg  
(Yeast)/ 0.02 mg (Baculovirus)/ 1 mg (E-Coli)/ 0.02 mg

(Mammalian-Cell)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg  
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-  
Cell)  
**NCBI Official Full Name**  
B3 domain-containing protein LOC\_Os02g10420

## NCBI Official Symbol

LOC107276623

## NCBI Official Synonym Symbols

OsJ\_05761; P0026H03.7; OJ1524\_D08.32

## NCBI Protein Information

B3 domain-containing protein LOC\_Os02g10420

## UniProt Gene Name

LOC\_Os02g10420

## UniProt Protein Name

B3 domain-containing protein LOC\_Os02g10420

## UniProt Primary Accession #

A3A463

## UniProt Secondary Accession #

Q6H744

## UniProt Related Accession #

A3A463

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**