

TYLCV C3 | Tomato yellow leaf curl virus coat protein C3

Cat PA00927

Size 150 μ l

Host

Rabbit

Clonality

Polyclonal

Confirmed reactivity

Tomato yellow leaf curl virus coat 16 kDa

Immunogen

Recombinant fragment of C3 protein as described in Gorovits et al. (2013).

Host

Rabbit

Clonality

Polyclonal

Purity

Serum

Format

Lyophilized

Reconstitution

For reconstitution add 150 μ l of sterile water

Storage

Store lyophilized/reconstituted at -20°C ; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Application

immunofluorescence (IF)

Recommended dilution

1 : 500 (IF)

Expected | apparent MW

16 kDa

Confirmed reactivity

Tomato yellow leaf curl virus coat 16 kDa

Predicted reactivity

Tomato yellow leaf curl virus coat 16 kDa

Not reactive in

No confirmed exceptions from predicted reactivity are currently known

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

TYLCV C3 | Tomato yellow leaf curl virus coat protein C3

Cat PA00927

Size 150 μ l

Additional information

This antibody is detecting recombinant TYLCV C3 protein

Description

TYLCV C3(Tomato Yellow Leaf Roll Virus coat protein C3) accumulates in tomato leaves during infection. The protein is detected in the phloem of the relevant cells. The TYLCV virus attracts tomato cultures worldwide and is transmitted by Whiteflies.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY