

Patatin-06, Recombinant Protein

Cat RP15383

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)

Full Product Name

Recombinant Solanum tuberosum Patatin-06

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

KLEEMVTVLS IDGGGIGKII PAIILEFLEG QLQEVNNDK ARLADYFDVI GGTSTGGLLT AMITTPNENN
RPFAAAKDIV PFYFEHGPHI FNYSGSIFGP RYDGKYLQV LQEKLGTRV HQALTEVAIS SFDIKTNKPV
IFTKSNLAKS PELDAKMYDI CYSTAAPIY FPPHHFVTHT SNGATYEFNL VGGVATVGD PALLSLSVAT
RLAQEDPAFS SIKSLDYKQM LLLSLGTGTN SEFDKTYTAE EAAKWGPLRW MLAIQMTNA ASSYMTDYYI
STVFQARHSQ NNYLRVQENA LTGTTTEMDD ASEANMELLV QVGETLLKVP VSKDSPETYE EALKRFAKLL
SNRKKLRANK ASY

Sequence Positions

24-386, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

42,614 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Patatin

NCBI Accession

Q2MY45.1

NCBI GI

122201880

NCBI Official Full Name

Patatin-06

UniProt Protein Name

Patatin-06

UniProt Primary Accession

Q2MY45

UniProt Comments

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Patatin-06, Recombinant Protein



Cat *RP15383*

Size *0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1*

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Probable lipolytic acyl hydrolase (LAP), an activity which is thought to be involved in the response of tubers to pathogens.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY