

Plant-Derived Horseradish Peroxidase Purified Lyophilized

Cat NP00045

Species

Horseradish Root

Purity

chromatographically purified

Format

lyophilized powder

Storage

2-8°C or -20°C

Concentration

500 Units/mg

Application

Plant-Derived Horseradish Peroxidase Purified lyophilized is useful as a hemoprotein that catalyzes oxidation by hydrogen peroxide of a number of substrates.

Description

Some examples of substrates include: ascorbate, ferrocyanide, cytochrome c, and the white form of many dyes. This applies to immune coupling.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY