

pGTQL1221YC

Cat PP00029

Vector type

Plant Expression ; BiFC assay

Standard format

Vector backbone

pEarleyGate-202

Modifications to backbone

YFP-C was fused after the Gateway cassette

Selectable markers

Basta

Bacterial Resistance(s)

Chloramphenicol and Kanamycin, 25 & 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DB3.1

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY