

pMpGWB329

Cat PP00122

Vector type

Plant Expression

Standard format

Vector backbone

pMpGWB300

Modifications to backbone

Coding sequence for tdTomato is cloned downstream of the Gateway cassette.

Selectable markers

Chlorsulfuron

Bacterial Resistance(s)

Chloramphenicol and Spectinomycin, 25 & 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DB3.1

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY