
Vector type

Plant Expression

Standard format

Vector backbone

pMpGWB300

Modifications to backbone

CaMV 35S promoter and coding sequence for tdTomato are cloned upstream and downstream of the Gateway cassette, respectively.

Selectable markers

Chlorsulfuron

Bacterial Resistance(s)

Chloramphenicol and Spectinomycin, 25 & 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DB3.1

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY