

---

## Vector type

Plant Expression

## Standard format

## Vector backbone

pMpGWB300

## Modifications to backbone

Coding sequence for Emerald Luc fused with PEST sequence [ELuc(PEST)] is cloned downstream of the Gateway cassette.

## Selectable markers

Chlorsulfuron

## Bacterial Resistance(s)

Chloramphenicol and Spectinomycin, 25 & 50 µg/mL

## Growth Temperature

37°C

## Growth Strain(s)

DB3.1

## Copy number

High Copy

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**