

pMpGWB332

Cat PP00125

Vector type

Plant Expression

Standard format

Vector backbone

pMpGWB300

Modifications to backbone

MpHSP17.8A1 promoter is cloned upstream of the Gateway cassette.

Selectable markers

Chlorsulfuron

Bacterial Resistance(s)

Chloramphenicol and Spectinomycin, 25 & 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DB3.1

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY