

pMpGWB400

Cat PP00130

Vector type

Plant Expression

Standard format

Vector backbone

pGWB400

Modifications to backbone

The proNOS:nptII:terNOS marker cassette was replaced with pro35Sx2:nptII:ter35S.

Selectable markers

Neomycin (select with G418)

Bacterial Resistance(s)

Spectinomycin, 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DH5α

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY