

pMpGWB414

Cat PP00144

Vector type

Plant Expression

Standard format

Vector backbone

pMpGWB400

Modifications to backbone

CaMV 35S promoter and glucocorticoid receptor (GR) coding sequence are cloned upstream and downstream of the Gateway cassette, respectively.

Selectable markers

Neomycin (select with G418)

Bacterial Resistance(s)

Chloramphenicol and Spectinomycin, 25 & 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DB3.1

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY