

pMpGWB415

Cat PP00145

Vector type

Plant Expression

Standard format

Vector backbone

pMpGWB400

Modifications to backbone

Nuclear localization signal-tagged Citrine (Citrine-NLS) coding sequence is cloned downstream of the Gateway cassette.

Selectable markers

Neomycin (select with G418)

Bacterial Resistance(s)

Chloramphenicol and Spectinomycin, 25 & 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DB3.1

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY