

# pMpGWB418

Cat PP00148

---

## Vector type

Plant Expression

## Standard format

## Vector backbone

pMpGWB400

## Modifications to backbone

MpEF1 $\alpha$  promoter and coding sequence for the modified EAR motif plant-specific repression domain showing strong repression activity (SRDX) are cloned upstream and downstream of the Gateway cassette, respectively.

## Selectable markers

Neomycin (select with G418)

## Bacterial Resistance(s)

Chloramphenicol and Spectinomycin, 25 & 50  $\mu$ g/mL

## Growth Temperature

37°C

## Growth Strain(s)

DB3.1

## Copy number

High Copy

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**