

pMpGWB420

Cat PP00150

Vector type

Plant Expression

Standard format

Vector backbone

pMpGWB400

Modifications to backbone

Coding sequence for a fusion polypeptide between the modified EAR motif plant-specific repression domain showing strong repression activity (SRDX) and glucocorticoid receptor (GR) is cloned downstream of the Gateway cassette.

Selectable markers

Neomycin (select with G418)

Bacterial Resistance(s)

Chloramphenicol and Spectinomycin, 25 & 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DB3.1

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY