

---

## Vector type

Plant Expression

## Standard format

## Total vector size (bp)

6580

## Vector backbone

pGreenII S BAR

## Modifications to backbone

35S CaMV promoter, Ocs 3' corrected terminator

## Selectable markers

Basta

## Bacterial Resistance(s)

Spectinomycin, 50 µg/mL

## Growth Temperature

37°C

## Growth Strain(s)

DH5α

## Copy number

High Copy

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**