

pRGEB32

Cat PP00415

Vector type

Plant Expression, CRISPR

Standard format

Vector backbone

pRGEB31

Modifications to backbone

Replace the 35S promoter with rice UBI promoter plus 5'-UTR.

Selectable markers

Hygromycin

Bacterial Resistance(s)

Kanamycin, 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DH5α

Copy number

High Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY